

# LANDMARK VARSITY

## Small Animal Anesthesia Cabinet System



### POWER, SAFETY & MOBILITY

The Landmark Varsity combines the best of both worlds – exceptional anesthesia delivery through our Landmark Table Top Anesthesia System-VTA-1100- and a convenient cabinet for storage.

Vetland Medical designed the Landmark Varsity Cabinet System to be a convenient, robust, easily maneuverable complete anesthesia system with easy attachment to O<sub>2</sub> E-tanks and our e2 Oxygen Concentrator.



11.25" x 20.25" workspace

built-in mount for E-Tank & Oxygen Concentrator (sold separately)



Accepts any vaporizer

Shown with (optional) calibrated for life Vetland D19 vaporizer (with 10 year warranty)

18"D x 14"W x 3.75"H drawers

18"D x 14"W x 5.75"H drawer

#### VAR-5500 Standard Features

- VTA-1100 anesthesia system with storage shelf & 3 large Secure-Lock drawers
- E-Tank mount with regulator
- Oxygen concentrator mount
- **Delux Safe™** pop-off valve with thumb occlusion
- Proven human-quality (FDA type) valves
- Heat & moisture conserving coaxial absorber
- Quick disconnect fresh gas outlet for non-rebreathing applications
- 0-4 LPM oxygen flow control with protective stop
- Accessory kit\*

#### What Makes the VAR-5500 Great?

- Neat, balanced, compact design
- Large Stainless steel workspace
- 6" casters; front 2 are locking
- 1030 cc absorber with less trapped volume, less chance of "channeling" & quicker anesthetic adjustments
- Low flow kit available
- Pull-out tray (optional)
- Back cover to protect all tubing connections
- Corrosion proof construction
- Capable of accepting any type of vaporizer
- **24 x 7 x 365 technical support**
- **5 Year Warranty**

VAR-5500 Landmark VTA-1100 & VAR-CART

VAR-CART Rolling Cabinet Only

Made in U.S.A. 

**\*Included Accessory Kit:** • Adult F Breathing Circuit • Modified Jackson Rees Non-Rebreathing Circuit • 1 Bag EcoSorb • 2-Liter Breathing Bag • 3-Liter Breathing Bag • Operator's Manual