

VETLAND V1407, V1410 & V1412

Veterinary Critical Care Patient Monitors



THE INDUSTRY'S BEST NIBP & SpO₂ TECHNOLOGY

The V1400-Series Monitors feature Sun Tech NIBP and Masimo SpO₂ technology - the recognized leaders in their respective fields.



V1400-Series Standard Features

- Veterinary specific algorithms
- Parameters Measured:
ECG, SpO₂, NIBP, Resp, 2-Temp, 2-IBP
- ECG-3 lead with 5 filtering modes & 4 sweep rates
- 13 arrhythmia analysis
- 7 days tabular & graphic trend data storage
- AC or long li-ion battery operation
- LAN & wireless LAN connection to central station
- SmartCare capable

What Makes the V1400-Series Monitors Great?

- Portable, compact design
- Intuitive 7", 10" or 12" high resolution display with maximum 10 waveforms
- Validated industry's best NIBP technology and SpO₂ technology from Sun Tech and Masimo
- Optional EtCO₂ module from Respirationics
- Optional touch screen
- Easy software upgradability
- PhaseIn's multi-gas mainstream (IRMA AX+), sidestream (ISA AX+, OR+)

Scan to
Watch the Video



WWW.VETLANDMEDICAL.COM

VETLAND V1407, 1410 & V1412

Veterinary Patient Monitors continued...

INNOVATIVE, VERSATILE & ACCURATE

The Vetland V1407, V1410 & V1412 Multi-Parameter Veterinary Patient Monitors for Critical Care are the most intuitive, versatile and accurate monitors on the veterinary market today.

The compact design and high resolution display make them easy to use and understand. ECG, SpO₂, NIBP, Resp, 2 Temp and 2IBP configurations come standard with these monitors.

Masimo SpO₂
Sun Tech NIBP
Respironics EtCO₂



V1407 OPTIONS:

V1407-BASE	Vetland V1407 [includes ECG, SpO ₂ , NIBP, Temp & RR]
V1407-OPT-TOUC	Touch Screen
V1407-OPT-CO2S	Sidestream
V1407-CO2M	Mainstream
V1407-OPT-PRIN	Printer
V-OPT-CART	Rolling Stand Mount
V1407-OPT-WALL	Wall Mount
V1407-OPT-TEMP	Temp

V1407-OPT-1IBP	Invasive Blood Pressure Monitoring
V1407-OPT-VGA	VGA Output
V1407-ESRT	Vetland V1407 STD-2 [ECG, SPO ₂ , RR & Temp only]
V1407-SNT	Vetland V1407 STD-3 [SpO ₂ , NIBP & Temp only]
V1407-SN	Vetland V1407 STD-4 [SpO ₂ & NIBP only]
V1407-C	Vetland V1407-E [CO ₂ only]
V1412-OPT-MAN	Operator's Manual




V1410 OPTIONS:

V1410-BASE	Vetland V1410 [includes ECG, SpO ₂ , NIBP, 2Temp, 2IBP & RR]
V1410-OPT-TOUC	Touch Screen
V1410-OPT-CO2S	Sidestream
V1410-CO2M	Mainstream
V1410-OPT-PRIN	Printer
V-OPT-CART	Rolling Stand Mount
V1410-OPT-WALL	Wall Mount
V1410-OPT-VGA	VGA Output
V1412-OPT-MAN	Operator's Manual

V1412 OPTIONS:

V1412-BASE	Vetland V1412 [includes ECG, SpO ₂ , NIBP, RR, 2Temp & 2IBP]
V1412-OPT-TOUC	Touch Screen
V1412-OPT-CO2S	Sidestream
V1412-OPT-CO2M	Mainstream
V1412-OPT-PRIN	Printer
V-OPT-CART	Rolling Stand Mount
V1412-OPT-WALL	Wall Mount
V1412-OPT-VGA	VGA Output
V1410/12-OPT-Man	Operator's Manual
V-OPT-MGAS	Multigas Option
VA-9900-9000-00	Battery for V1410/V1412



WWW.VETLANDMEDICAL.COM