

VETLAND V1407, V1410 & V1412

Veterinary Anesthesia Critical Care Patient Monitors



THE INDUSTRY'S BEST NIBP & SpO₂ TECHNOLOGY

The V1400-Series Monitors feature Sun Tech NIBP and Masimo SpO₂ technology - the recognized leaders in their respective fields.



V1400-Series Standard Features

- Veterinary specific algorithms
- Color touchscreen display
- Parameters Measured:
ECG, SpO₂, NIBP, Resp, *2Temp, *2IBP
- ECG-3 lead with 5 filtering modes & 4 sweep rates
- 13 arrhythmia analysis
- 7 days tabular & graphic trend data storage
- AC or long li-ion battery operation
- LAN & wireless LAN connection to central station
- AAHA compliant data importability to any practice management software via SD card or network connection

What Makes the V1400-Series Monitors Great?

- Portable, compact design
- Color touchscreen
- Intuitive 7", 10" or 12" high resolution display with maximum 10 waveforms
- Validated industry's best NIBP technology and SpO₂ technology from Sun Tech and Masimo
- Optional EtCO₂ module from Respirationics
- Easy software upgradability
- Optional remote viewing from anywhere in the world via Smart Care System
- Optional PhaseIn's multi-gas mainstream (IRMA AX+, sidestream (ISA AX+, OR+)

* 1410 & 1412 models only



WWW.VETLANDMEDICAL.COM

SCAN
TO
LEARN
MORE



VETLAND V1407, V1410 & V1412

Veterinary Anesthesia Critical Care Patient Monitors continued...

INNOVATIVE, VERSATILE & ACCURATE

The Vetland V1407, V1410 & V1412 Multi-Parameter Veterinary Patient Monitors for Critical Care are the most intuitive, versatile and accurate monitors on the veterinary market today.

The compact design and high resolution display make them easy to use and understand.

ECG, SpO₂, NIBP, Resp, *2Temp and *2IBP configurations come standard with these monitors.

* 1410 & 1412 models only

Masimo
SpO₂ SunTech
NIBP Respiration
EtCO₂



V1407 OPTIONS:

V1407-BASE	Vetland V1407 [includes ECG, SpO ₂ , NIBP, Temp & RR]
V-OPT-CO2S	Sidestream
V-OPT-CO2M	Mainstream
V-OPT-PRIN	Printer
V-OPT-CART	Rolling Stand Mount
V-OPT-WALL	Wall Mount
V1407-OPT-IBP	Invasive Blood Pressure Monitoring
V-OPT-VGA	VGA Output
V1407-ESRT	Vetland V1407 STD-2 [ECG, SPO ₂ , RR & Temp only]
V1407-SNT	Vetland V1407 STD-3 [SpO ₂ , NIBP & Temp only]
V1407-SN	Vetland V1407 STD-4 [SpO ₂ & NIBP only]
V1407-C	Vetland V1407-E [CO ₂ only]
V1407-OPT-MAN	Operator's Manual

V1410 OPTIONS:

V1410-BASE	Vetland V1410 [includes ECG, SpO ₂ , NIBP, 2Temp, 2IBP & RR]
V-OPT-CO2S	Sidestream
V-CO2M	Mainstream
V-OPT-PRIN	Printer
V-OPT-CART	Rolling Stand Mount
V-OPT-WALL	Wall Mount
V-OPT-VGA	VGA Output
VA-9900-9000-00	Battery for V1410/V1412
V1410/12-OPT-MAN	Operator's Manual

V1412 OPTIONS

V1412-BASE	Vetland V1412 [includes ECG, SpO ₂ , NIBP, RR, 2Temp & 2IBP]
V-OPT-CO2S	Sidestream
V-OPT-CO2M	Mainstream
V-OPT-PRIN	Printer
V-OPT-CART	Rolling Stand Mount
V-OPT-WALL	Wall Mount
V-OPT-VGA	VGA Output
VA-9900-9000-00	Battery for V1410/V1412
V-OPT-MGAS	Multigas Option
V1410/12-OPT-MAN	Operator's Manual



WWW.VETLANDMEDICAL.COM